




www.biocare.com.cn

BP-12C

Handheld Pulse
Oximeter



Features:

- Elegant and handheld design
- Color TFT LCD display of SpO₂, PR, Pulse bar, PI and Plethysmogram
- Real-time review SpO₂ and PR value
- Up to 384-hour of oximeter data storage (Record data once every 1/2/4/8 seconds)
- Oximeter data manager software is available
- Audible & visible alarm function
- Low battery voltage indication
- Power off automatically when no signal
- More than 40 hours battery life
- AC/DC adapter (optional)
- Data transmission to PC for view and analysis (optional)



Specifications

Model		BP-12C
SpO ₂	Measurement range	35%-100%
	Accuracy	≤ 3% (SpO ₂ range from 70% ~ 100%)
	SpO ₂ low over-limit setting range	50%~99%
	Default low over-limit	90%(adult), 95%(pediatric)
PR	Measurement range	30 bpm~240 bpm
	Accuracy	± 2% or ± 2 bpm (whichever is greater)
	PR over-limit setting range	25 bpm~250 bpm
	Default over-limit	120 bpm (high), 50 bpm (low)
Perfusion Index	Measurement range	0.2%~20%
Power supply	Battery type	3 x LR6(AA)alkaline batteries or Ni-MH rechargeable batteries(optional)
	Supply voltage	D.C. 4.5V (3 AA alkaline batteries),or D.C. 3.6V(Lithium rechargeable battery)
	Operating current	≤ 180mA
	Continuous working time	≤ 40 hours
Waterproof		IP 22
Data Record		Record data once every 1/2/4/8 seconds, up to 384-hour records
Dimensions		145 mm (L) × 74 mm (W) × 29 mm (H)
Net Weight		210g (including batteries)

Shenzhen Biocare Bio-Medical Equipment Co., Ltd.

Add: #16-1, Jinhui Road, Jinsha Community, Kengzi Sub-District, Pingshan New District, 518122 Shenzhen, P.R. CHINA

Tel: +86-755-2796 0888/3300 5899

Fax: +86-755-2796 0643

E-mail: sales@biocare.com.cn

Website: www.biocare.com.cn

© 2019 BIO CARE all rights reserved. Specifications subject to changes without prior notice.

ENG-CATA-V2.3-20190815

