



NEW POCKET SONOVIT

New Pocket Sonovit is the latest professional device for ultrasound therapy. It is easy to use and allows many applications. Both probe and device, in fact, are able to supply vibrations for the execution of various therapeutic applications. New Pocket Sonovit contains pre-set programmes for the more utilised applications; functions by power mains and with an internal rechargeable battery. New Pocket Sonovit is equipped with a handy bag case for transport and home therapy.



Programmes / Applications

Bone Atrophy	Muscle Tears
Bone Callus	Contractures
Osteitis and Periostitis	Bruises
Delayed Ossification	Tenosynovitis
Arthritis / Arthrosis	Epicondylitis
Joint Stiffness	Periarthritis
Sprains	Bursitis
Lumbago	Tendonitis
Myalgia	Cicatrisation
Spasms	

Technical Specifications

Power Supply	rechargeable battery / power mains
No. Output Channels	1
No. of Programmes	20 pre-set at 3MHz + 10 free at 1/3MHz
Frequency (Hz)	1/3 MHz
Emission	continuous, modulated (1/2, 1/4, 1/8, 1/10)
Power	3 W/cm ²
Dimensions and weight	18*11*4cm, 400 gr only device, 1,5 kg with accessories

Standard Kit

Bag
 Ultrasound probe Ø 45 mm
 Rechargeable battery
 Power supply
 User's manual with photos of the applications